

Product Datasheet

MRS-1706 (orb1218249)

Catalog Number	orb1218249
Description	MRS-1706 has Ki values of 1.39, 112, 157, and 230 nM for human A2B, A2A, A1, and A3 receptors respectively. MRS-1706 is an effective and selective adenosine A2B receptor inverse agonist.
Target	A1 A2A A2B A3
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{27}H_{29}N_5O_5$
Purity	> 98% (HPLC)
MW	503.6
CAS Number	264622-53-9
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)