

## Product Datasheet

### Solabegron (orb1218279)

<b>Catalog Number</b>	orb1218279
<b>Category</b>	Small Molecules
<b>Description</b>	Solabegron is a selective agonist of $\beta$ 3-adrenergic receptor, stimulating cAMP accumulation in Chinese hamster ovary cells expressing the human $\beta$ 3-AR(EC50 : 22 nM). It is used for treatment of overactive bladder and irritable bowel syndrome.
<b>Target</b>	$\beta$ 3-adrenergic receptor
<b>Purity</b>	>98% (HPLC)
<b>MW</b>	410.9
<b>CAS Number</b>	252920-94-8
<b>Formula</b>	$C_{23}H_{23}ClN_2O_3$
<b>Storage</b>	Storage temperature: -20°C. Stability: $\geq$ 2 years
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)