

Product Datasheet

WUN29654 (orb1218318)

Catalog Number	orb1218318
Category	Small Molecules
Description	WUN29654 is a Thalidomide analog and CRBN modulator, has an IC ₅₀ of 3.5 μM and a K _i of 0.98 μM.
Target	Others
Purity	>98% (HPLC)
MW	324.29
Solubility (25°C)	DMSO:10 mM
CAS Number	2407829-65-4
Formula	C ₁₇ H ₁₂ N ₂ O ₅
SMILES	<chem>O=C(C1=C2C(C=C(O)C=C23)=CC=C1)N(C4CCC(NC4=O)=O)C3=O</chem>
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558