

Product Datasheet

Ac4ManNAz (orb1220165)

Catalog Number	orb1220165
Category	Small Molecules
Description	Ac4ManNAz is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs.
Target	PEGs
Purity	>98% (HPLC)
MW	430.4
CAS Number	361154-30-5
Formula	C ₁₆ H ₂₂ N ₄ O ₁₀
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)