

Product Datasheet

Protease-Activated Receptor-4 (orb1220396)

Catalog Number	orb1220396
Description	Protease-Activated Receptor-4 is the proteinase-activated receptor-4 (PAR4) agonist.
Target	proteinase-activated receptor-4 (PAR-4)
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{33}H_{46}N_8O_7$
Purity	> 98% (HPLC)
MW	666.8
CAS Number	245443-52-1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)