

Product Datasheet

Yoda 1 (orb1220408)

Catalog Number	orb1220408
Description	Yoda1 is a Piezo1 ion channel agonist (EC50s = 17.1 and 26.6 μ M for murine and human Piezo1-transfected HEK293T cells respectively).
Target	Piezo1
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only
Formula	$C_{13}H_8Cl_2N_4S_2$
Purity	> 98% (HPLC)
MW	355.3
CAS Number	448947-81-7
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)