

Product Datasheet

F-15599 (orb1220487)

Catalog Number	orb1220487
Description	F-15599 is a highly selective G-protein biased 5-HT1A receptor agonist(Ki : 3.4 nM).
Target	5-HT1A
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{19}H_{21}ClF_2N_4O$
Purity	> 98% (HPLC)
MW	394.9
CAS Number	635323-95-4
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)