

Product Datasheet

GSK-4716 (orb1220987)

Catalog Number	orb1220987
Description	GSK4716 is an ERRβ/γ agonist (IC ₅₀ = 2 μM) that can activate the receptor with similar potency as the protein ligand.
Target	ERRβ/γ
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	C ₁₇ H ₁₈ N ₂ O ₂
Purity	> 98% (HPLC)
MW	282.3
CAS Number	101574-65-6
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)