

Product Datasheet

Treprostinil (orb1221185)

Catalog Number	orb1221185
Description	Treprostinil is a potent DP1 IP and EP2 agonist (EC50: 0.6/1.9/6.2 nM).
Target	DP1 IP EP2
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{23}H_{34}O_5$
Purity	> 98% (HPLC)
MW	390.5
CAS Number	81846-19-7
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)