

Product Datasheet

Clonidine (orb1221189)

Catalog Number	orb1221189
Description	Clonidine is a centrally active alpha-adrenergic agonist used predominantly as an antihypertensive agent usually in combination with other agents. Despite wide-scale use for many years clonidine has not been linked definitively to either serum aminotransferase elevations or clinically apparent liver injury.
Target	α 2 adrenergic receptor
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only
Formula	$C_9H_9Cl_2N_3$
Purity	> 98% (HPLC)
MW	230.1
CAS Number	4205-90-7
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)