

Product Datasheet

Isoetharine mesylate salt (orb1221213)

Catalog Number	orb1221213
Description	Isoetharine Mesylate Salt is a β -adrenergic receptor agonist. It also is pregnane X receptor (PXR) activator capable of upregulating CYP450 expression.
Target	β -adrenergic receptor PXR
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only
Formula	$C_{14}H_{25}NO_6S$
Purity	> 98% (HPLC)
MW	335.4
CAS Number	7279-75-6
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)