

Product Datasheet

Uridine 5'-triphosphate tris salt (orb1221375)

Catalog Number	orb1221375
Description	Uridine 5'-triphosphate tris salt is an agonist to the P2Y receptor and has been used to study concentrations of intracellular Ca ²⁺ as well as a phosphate donor for in vitro kinase assays.
Target	Others
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{13}H_{26}N_3O_{18}P_3$
Purity	> 98% (HPLC)
MW	605.3
CAS Number	108321-53-5
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)