

Product Datasheet

Quinolinic acid (orb1221680)

Catalog Number	orb1221680
Category	Small Molecules
Description	Quinolinic acid is an endogenous N-methyl-D-aspartate (NMDA) receptor agonist synthesized from L-tryptophan via the kynurenine pathway.
Target	NMDA receptor
Purity	>98% (HPLC)
MW	173.2
CAS Number	89-00-9
Formula	$C_{10}H_7NO_2$
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)