

Product Datasheet

Succinylcholine Chloride Dihydrate (orb1221968)

Catalog Number	orb1221968
Description	Succinylcholine Chloride Dihydrate is a nicotinic AChR agonist and also acts as a depolarizing neuromuscular blocker.
Target	nicotinic AChR
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{14}H_{30}N_2O_4 \cdot 2Cl \cdot 2H_2O$
Purity	> 98% (HPLC)
MW	397.3
CAS Number	6101-15-1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)