

Product Datasheet

TLR7-agonist-1 (orb1222628)

Catalog Number	orb1222628
Description	TLR7-agonist-1 is a specific and effective Toll-like Receptor 7 (TLR7) agonist (LEC=0.4 μ M).
Target	TLR7
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only
Formula	$C_{17}H_{16}N_6O_2$
Purity	> 98% (HPLC)
MW	336.4
CAS Number	1642857-69-9
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)