

Product Datasheet

4-IBP (orb1222643)

Catalog Number	orb1222643
Description	4-IBP is a selective σ_1 agonist with a high level of affinity for the σ_1 receptor (K_i = 1.7 nM) and a moderate affinity for the σ_2 receptor (K_i = 25.2 nM) .
Target	$\sigma_1 \sigma_2$
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only
Formula	$C_{19}H_{21}IN_2O$
Purity	> 98% (HPLC)
MW	420.3
CAS Number	155798-08-6
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)