

Product Datasheet

Adenosine 5'-triphosphate disodium salt (orb1223480)

Catalog Number	orb1223480
Description	ATP Disodium salt is a P2 purinoceptor agonist.
Target	Others
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{10}H_{14}N_5O_{13}P_3 \cdot 2Na$
Purity	> 98%(HPLC)
MW	551.1
CAS Number	987-65-5
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)