

Product Datasheet

Pimavanserin (orb1224181)

Catalog Number	orb1224181
Category	Small Molecules
Description	Pimavanserin(ACP-103) is a potent and selective 5-HT2A receptor inverse agonist with mean pIC50 of with 8.7 in the cell-based functional assay.
Target	5-HT2A
Purity	>98%(HPLC)
MW	427.6
CAS Number	706779-91-1
Formula	$C_{25}H_{34}FN_3O_2$
Storage	Storage temperature: -20°C. Stability: \geq 2 years
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)