

Product Datasheet

Brimonidine Tartrate (orb1224186)

Catalog Number	orb1224186
Description	Brimonidine Tartrate is a highly selective α -adrenergic receptor agonist with EC ₅₀ of 0.45 nM for the α 2A adrenoreceptor, and used to treat open-angle glaucoma or ocular hypertension.
Target	α 2A-adrenergic receptor
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only
Formula	C ₁₅ H ₁₆ BrN ₅ O ₆
Purity	> 98%(HPLC)
MW	442.2
CAS Number	70359-46-5
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)