

Product Datasheet

Pergolide mesylate salt (orb1224311)

Catalog Number	orb1224311
Description	PERGOLIDE MESYLATE is a long-acting dopamine agonist which has been used to treat parkinson disease and hyperprolactinemiabut withdrawn from some markets due to potential for heart valve diseases.
Target	5-HT receptor Adrenergic receptor Dopamine receptor
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{20}H_{30}N_2O_3S_2$
Purity	> 98%(HPLC)
MW	410.6
CAS Number	66104-23-2
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)