

Product Datasheet

Cyproterone acetate (orb1225070)

Catalog Number	orb1225070
Description	Cyproterone acetate is an androgen receptor (AR) antagonist with IC50 of 7.1 nM, as well as a weak progesterone receptor agonist with weak pro-gestational and glucocorticoid activity.
Target	Androgen Receptor
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{24}H_{29}ClO_4$
Purity	> 98%(HPLC)
MW	416.9
CAS Number	427-51-0
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)