

Product Datasheet

Pardoprunox (orb1225462)

Catalog Number	orb1225462
Description	Pardoprunox(SLV-308) is a novel partial dopamine D2 and D3 receptor agonist and serotonin 5-HT1A receptor agonist.
Target	5-HT1A α 2-adrenergic receptor D2 D3
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only
Formula	$C_{12}H_{15}N_3O_2$
Purity	> 98%(HPLC)
MW	233.3
CAS Number	269718-84-5
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)