

Product Datasheet

WHI-P180 (orb1225684)

Catalog Number	orb1225684
Category	Small Molecules
Description	WHI-P180 is a potent EGFR and Cdk2 inhibitors with IC50 of 4.0 and 1.0 uM, respectively.
Target	CDK2 EGFR
Purity	>98%(HPLC)
MW	297.3
CAS Number	211555-08-7
Formula	$C_{16}H_{15}N_3O_3$
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)