

Product Datasheet

IL-23R Recombinant Protein (orb1228390)

| | |
|--------------------------|---|
| Catalog Number | orb1228390 |
| Description | IL-23R Recombinant Protein |
| Target | IL23r |
| Preservatives | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O. |
| Form/Appearance | Lyophilized |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | C-Fc tag |
| Note | For research use only |
| Application notes | This recombinant protein can be used for biological assays. For research use only. |
| Protein Sequence | Gly24-Asp372 |
| Purity | Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test. |
| Source | Human Cells |
| MW | 67.3 kD |
| Uniprot ID | Q8K4B4 |
| NCBI | Q8K4B4 |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)