

Product Datasheet

NCOA4 (IN) Peptide (orb1231635)

Catalog Number	orb1231635
Description	NCOA4 (IN) Peptide
Tested Applications	Blocking
Target	NCOA4
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store NCOA4 (IN) Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	NCOA4 (IN) Peptide is used for blocking the activity of the NCOA4 (IN) Antibody
Uniprot ID	Q13772
NCBI	Q13772
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)