

Product Datasheet

USP25 Peptide (orb1231923)

Catalog Number	orb1231923
Category	Proteins
Description	USP25 Peptide
Target	USP25
Form/Appearance	Liquid
Concentration	200 µg/mL
Buffer/Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
UniProt ID	Q9UHP3
Tested applications	Blocking
Application notes	The peptide is used for blocking the activity of USP25 Antibody. It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Storage	Store USP25 Peptide at -20°C, stable for one year.
Note	For research use only
NCBI	NP_037528
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Western blot analysis of USP25 in mouse brain tissue lysate with USP25 antibody at 1 $\mu\text{g}/\text{ml}$ in (A) the absence and (B) the presence of blocking peptide.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)