

Product Datasheet

ABCA7 Peptide (orb1232218)

Catalog Number	orb1232218
Category	Proteins
Description	ABCA7 Peptide
Target	ABCA7
Form/Appearance	Liquid
Concentration	200 µg/mL
Buffer/Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
UniProt ID	Q8IZY2
Tested applications	Blocking
Application notes	ABCA7 peptide is used for blocking the activity of ABCA7 antibody.
Storage	Store ABCA7 peptide at -20°C, stable for one year.
Note	For research use only
NCBI	NP_061985
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)