

Product Datasheet

E2-EPF Peptide (orb1234209)

Catalog Number	orb1234209
Category	Proteins
Description	E2-EPF Peptide
Target	UBE2S
Form/Appearance	Lyophilized
Buffer/Preservatives	E2-EPF peptide is supplied as 100 µg of dried peptide. When the peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL.
UniProt ID	Q16763
Tested applications	Blocking
Application notes	E2-EPF Peptide is used for blocking the activity of E2-EPF antibody. Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been pre-absorbed with the provided peptide.
Storage	Store at -20°C. When reconstituted minimize freezing and thawing.
Note	For research use only
NCBI	NP_055316.2
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)