

Product Datasheet

Transcription factor E2F4 Peptide (orb1235055)

Catalog Number	orb1235055
Category	Proteins
Description	Transcription factor E2F4 Peptide
Target	Transcription factor E2F4
Form/Appearance	Lyophilized
Buffer/Preservatives	Transcription factor E2F4 peptide is supplied as 100 µg of dried peptide. When the peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL.
UniProt ID	Q16254
Tested applications	Blocking
Application notes	Transcription factor E2F4 peptide is used for blocking the activity of Transcription factor E2F4 antibody. Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been pre-absorbed with the provided peptide.
Storage	Store at -20°C. When reconstituted minimize freezing and thawing.
Note	For research use only
NCBI	NP_001941.2
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)