

Product Datasheet

PCBP1 Peptide (orb1235763)

Catalog Number	orb1235763
Description	PCBP1 Peptide
Tested Applications	Blocking
Target	PCBP1
Preservatives	PCBP1 peptide is supplied as 100 µg of dried peptide. When the peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL.
Form/Appearance	Lyophilized
Storage	Store at -20°C. When reconstituted minimize freezing and thawing.
Note	For research use only
Application notes	PCBP1 Peptide is used for blocking the activity of PCBP1 antibody. Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been pre-absorbed with the provided peptide.
Uniprot ID	Q15365
NCBI	NP_006187.2
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)