

Product Datasheet

Donkey Anti-Goat IgG (H&L) Antibody (Texas Red) (orb1236949)

Description	Donkey Anti-Goat IgG (H&L) Antibody (Texas Red)
Species/Host	Donkey
Reactivity	Goat
Conjugation	Texas Red
Tested Applications	ELISA, IHC, WB
Immunogen	Goat IgG - whole molecule
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, 0.01% (w/v) Sodium Azide; pH 7.2; 10 mg/mL Bovine Serum Albumin
Form/Appearance	Lyophilized
Storage	Alkaline Phosphatase Conjugated Affinity Purified anti-Rabbit IgG can be stored vial at 4 °°C before opening. DO NOT FREEZE. This product is stable at 4 °°C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity. Expiration date is one (1) year from date of opening product.
Note	For research use only
Application notes	Texas Red Conjugated Affinity Purified anti-Goat IgG is suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency.
Isotype	IgG
Antibody Type	Secondary Antibody
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)