

Product Datasheet

Mouse Anti-Goat IgG Antibody [1G12H5] (HRP) (orb1236981)

Description	Mouse Anti-Goat IgG Antibody [1G12H5] (HRP)
Species/Host	Mouse
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Immunogen	Goat IgG
Preservatives	PBS pH 7.2, 1% rAlbumin, 0.02% Thimerosal
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 2-8°C. The product is stable for about 4-6 weeks at 4°C as an undiluted liquid. For extended storage, aliquot and freeze the product at -70°C or below. Avoid repeated freezing and thawing. Or alternatively, add an equal volume of glycerol for a final concentration of 50% and store at -20°C.
Note	For research use only
Application notes	Immunoblotting: 1/5,000 - 1/50,000, ELISA: 1/5,000 - 1/50,000, Immunohistochemistry: 1/100 - 1/1,000
Isotype	IgG2
Clone Number	[1G12H5]
Antibody Type	Secondary Antibody
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)