

Product Datasheet

Control Normal Rabbit IgG (PE) Antibody (orb1237114)

Catalog Number	orb1237114
Category	Antibodies
Description	Control Normal Rabbit IgG (PE) Antibody
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Form/Appearance	Liquid
Concentration	batch dependent
Buffer/Preservatives	Supplied in PBS/sodium azide and a stabilizing agent
Tested applications	FC, FLISA
Application notes	Normal Rabbit IgG (PE) can be used as control in flow cytometry assays.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)