

Product Datasheet

GAG Antibody (HRP) (orb1238362)

Catalog Number	orb1238362
Description	GAG Antibody (HRP)
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Mouse monoclonal HIV-1 p24 antibody was raised against a recombinant full-length HIV-1 p24 protein.
Target	GAG
Preservatives	HIV-1 p24 antibody is supplied in PBS containing 1% rAlbumin and 0.02% thimerosol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	orb1238362 can be used for detection of p24 by Western blot or ELISA at 0.2 – 0.5µg/mL.
Isotype	IgG1
Clonality	Monoclonal
Clone Number	8G9
Antibody Type	Primary Antibody


Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

MW	Predicted: 24, 41, 55 kDa
Uniprot ID	P04591
NCBI	AAB50258
Expiration Date	12 months from date of receipt.



Western blot analysis of 20 ng of recombinant HIV-1 p24 protein with orb1238362 at 0.2 µg/mL.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)