

Product Datasheet

SARS-CoV-2 (COVID-19) ORF9c Antibody (orb1239056)

Catalog Number	orb1239056
Category	Antibodies
Description	SARS-CoV-2 (COVID-19) ORF9c Antibody
Target	ORF9c
Clonality	Polyclonal
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Virus
Form/Appearance	Liquid
Concentration	1 mg/mL
Buffer/Preservatives	SARS-CoV-2 (COVID-19) ORF9c antibody is supplied in PBS containing 0.02% sodium azide.
Purification	SARS-CoV-2 (COVID-19) ORF9c Antibody is affinity chromatography purified via peptide column.
Immunogen	Anti-SARS-CoV-2 (COVID-19) ORF9c antibody (orb1239056) was raised against a peptide corresponding to 16 amino acids near the amino-terminus of SARS-CoV-2 (COVID-19) ORF9c protein. The immunogen is located within the first 50 amino acids of the SARS-CoV-2 (COVID-19) ORF9c protein.
UniProt ID	P0DTD3
MW	8kD


Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)


Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Tested applications	ELISA, WB
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	P0DTD3



Western Blot Validation with SARS-CoV-2 (COVID-19) ORF9c Overexpressed 293 Cells. Loading: 10 µg per lane of 293 cell lysate from WT and SARS-CoV-2 (COVID-19) ORF9c transfected cells. Antibodies: SARS-CoV-2 (COVID-19) ORF9c orb1239056, 1 µg/mL, 1 h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.




Western Blot Validation with SARS-CoV-2 (COVID-19) ORF9c Overexpressed 293 Cells. Loading: 10 µg per lane of 293 cell lysate from WT and SARS-CoV-2 (COVID-19) ORF9c transfected cells. Antibodies: SARS-CoV-2 (COVID-19) ORF9c, orb1239056, 1 h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution. Lane 1: 0.25 µg/mL and. Lane 2: 0.5 µg/mL.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



ELISA Validation. Antibodies: SARS-CoV-2 (COVID-19) ORF9c Antibody, orb1239056. A direct ELISA was performed using SARS-CoV-2 ORF9c immunogen peptide (orb1231666) as coating antigen and the anti-SARS-CoV-2 (COVID-19) ORF9c antibody as the capture antibody. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:20000 dilution. Detection range is from 0.3 ng/mL to 1000 ng/mL.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)