

Product Datasheet

CORO7 Antibody (orb1239960)

Catalog Number	orb1239960
Category	Antibodies
Description	CORO7 Antibody
Target	CORO7
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine
Form/Appearance	Liquid
Concentration	1 mg/mL
Buffer/Preservatives	Coronin 7 Antibody is supplied in PBS containing 0.02% sodium azide.
Purification	Coronin 7 Antibody is affinity chromatography purified via peptide column.
Immunogen	Coronin 7 antibody was raised against an 17 amino acid synthetic peptide near the carboxy terminus of human Coronin 7.The immunogen is located within the last 50 amino acids of Coronin 7.
UniProt ID	P57737
Tested applications	ELISA, IF, IHC-P, WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Application notes

Coronin 7 antibody can be used for detection of Coronin 7 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2.5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Specificity

Multiple isoforms of Coronin 7 are known to exist

Antibody Type

Primary Antibody

Modifications

None

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only

NCBI

NP_078811

Expiration Date

12 months from date of receipt.



Western blot analysis of Coronin 7 in rat lung tissue lysate with Coronin 7 antibody at (A) 1 and (B) 2 µg/mL.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Immunohistochemistry of Coronin 7 in rat lung tissue with Coronin 7 antibody at 2.5 µg/mL.



Immunofluorescence of Coronin 7 in rat lung tissue with Coronin 7 antibody at 20 µg/mL.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)