

Product Datasheet

CUL9 Antibody (orb1249730)

Catalog Number	orb1249730
Category	Antibodies
Description	CUL9 Antibody
Target	CUL9
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	85 mg/ml
Buffer/Preservatives	Sterile filtered antiserum, 0.01% sodium azide.
Purification	Delipidation and Defibrination
Immunogen	Cullin 9 antibody was raised against amino acids 2503 - 2517 of Cullin 9 (Human).
UniProt ID	Q8IWT3
Tested applications	ELISA, IHC, IP, WB
Application notes	Cullin 9 antibody can be used in ELISA starting at 1:000 - 1:1000, and immunohistochemistry starting at 10 µg/mL.
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Modifications	None
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	Q8IWT3
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of Cullin 9 in colon tissue using Cullin 9 Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)