

Product Datasheet

LAMP2 Antibody (FITC) (orb1250585)

Catalog Number	orb1250585
Category	Antibodies
Description	LAMP2 Antibody (FITC)
Target	LAMP2
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	FITC
Reactivity	Human, Monkey, Primate
Form/Appearance	Liquid
Concentration	batch dependent
Buffer/Preservatives	Supplied in PBS/NaN ₃
Immunogen	Human adherent peripheral blood cells
UniProt ID	P13473
Tested applications	FC
Application notes	CD107b Antibody [H4B4] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and western blotting assays.
Specificity	CD107b

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Primary Antibody
Clone Number	[H4B4]
Modifications	None
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	P13473
Expiration Date	12 months from date of receipt.



Human T lymphocyte cell line Jurkat was intracellularly stained with Mouse Anti-Human CD107b-FITC (Cat. No. orb1250585).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)