

Product Datasheet

S-Tag Antibody [SBSTAGa] (FITC) (orb1251045)

Catalog Number	orb1251045
Category	Antibodies
Description	S-Tag Antibody [SBSTAGa] (FITC)
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	FITC
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Buffer/Preservatives	Supplied in PBS/NaN ₃
Immunogen	Recombinant S-tag protein
Tested applications	FC, IHC-Fr, IP
Application notes	S-Tag Antibody [SBSTAGa] for use in western blotting and other fluorescent assays.
Specificity	S-Tag
Antibody Type	Primary Antibody
Clone Number	[SBSTAGa]
Modifications	None

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Store vial at 2-8°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)