

Product Datasheet

Lamp1 Antibody (APC) (orb1251595)

Catalog Number	orb1251595
Category	Antibodies
Description	Lamp1 Antibody (APC)
Target	Lamp1
Clonality	Monoclonal
Species/Host	Rat
Isotype	IgG2a
Conjugation	APC
Reactivity	Mouse
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Buffer/Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Immunogen	Plasma membrane fraction of mouse embryo NIH/3T3 cell line
UniProt ID	P11438
Tested applications	Blocking, Cell Separation, FC, ICC, IHC-Fr, IP
Application notes	CD107a Antibody [1D4B] for use in flow cytometry assays.
Specificity	CD107a
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	[1D4B]
Modifications	None
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	P11438
Expiration Date	12 months from date of receipt.



Mouse cells were intracellularly stained with Rat Anti-Mouse CD107a-APC (Cat. No. orb1251595).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558