

Product Datasheet

CD16/32 Antibody [93] (PE-Cy5.5) (orb1251735)

Catalog Number	orb1251735
Description	CD16/32 Antibody [93] (PE-Cy5.5)
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Adhesion, Blocking, Depletion, ELISA, FC, ICC, IHC-Fr, IP
Immunogen	VH81X Tg B cells and T220 stromal cell line
Target	Fcgr2b; Fcgr3
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	CD16/32 Antibody [93] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Clone Number	[93]
Antibody Type	Primary Antibody
Uniprot ID	P08101, P08508

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

NCBI**P08508, P08101****Expiration Date**

12 months from date of receipt.



BALB/c mouse splenocytes were stained with Rat Anti-Mouse CD16/32-PE/CY5.5 (Cat. No. orb1251735) and Rat Anti-Mouse CD19-PE.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)