

## Product Datasheet

### PAPSS2 Antibody (orb1262760)


|                            |                                                                                                                                                                          |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | PAPSS2 Antibody                                                                                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                             |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                |
| <b>Immunogen</b>           | This PAPSS2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 544-574 amino acids from the C-terminal region of human PAPSS2. |
| <b>Target</b>              | PAPSS2                                                                                                                                                                   |
| <b>Preservatives</b>       | Supplied in PBS with 0.09% (W/V) sodium azide.                                                                                                                           |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                   |
| <b>Concentration</b>       | batch dependent                                                                                                                                                          |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.                                  |
| <b>Note</b>                | For research use only                                                                                                                                                    |
| <b>Application notes</b>   | For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100                                                                                              |
| <b>Isotype</b>             | Rabbit Ig                                                                                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                               |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                                                         |
| <b>MW</b>                  | 70 kDa                                                                                                                                                                   |

#### Biorbyt Ltd.


7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Uniprot ID** **O95340****NCBI** **O95340****Expiration Date** 12 months from date of receipt.

Antibody is used in Western blot to detect PAPSS2 in Jurkat cell lysate.



Formalin-fixed and paraffin-embedded human cancer tissue reacted with the primary antibody, which was peroxidase-conjugated to the secondary antibody, followed by AEC staining. BC = breast carcinoma; HC = hepatocarcinoma.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)