

## Product Datasheet

### Ntrk3 Antibody (orb1268172)

|                             |                                                                                                                                                                              |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1268172                                                                                                                                                                   |
| <b>Category</b>             | Antibodies                                                                                                                                                                   |
| <b>Description</b>          | Ntrk3 Antibody                                                                                                                                                               |
| <b>Target</b>               | Ntrk3                                                                                                                                                                        |
| <b>Clonality</b>            | Polyclonal                                                                                                                                                                   |
| <b>Species/Host</b>         | Rabbit                                                                                                                                                                       |
| <b>Isotype</b>              | Rabbit Ig                                                                                                                                                                    |
| <b>Conjugation</b>          | Unconjugated                                                                                                                                                                 |
| <b>Reactivity</b>           | Mouse                                                                                                                                                                        |
| <b>Predicted Reactivity</b> | Rat                                                                                                                                                                          |
| <b>Form/Appearance</b>      | Liquid                                                                                                                                                                       |
| <b>Concentration</b>        | batch dependent                                                                                                                                                              |
| <b>Buffer/Preservatives</b> | Supplied in PBS with 0.09% (W/V) sodium azide.                                                                                                                               |
| <b>Purification</b>         | This antibody is purified through a protein A column, followed by peptide affinity purification.                                                                             |
| <b>Immunogen</b>            | This Mouse Ntrk3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 161-188 amino acids from the N-terminal region of mouse Ntrk3. |
| <b>UniProt ID</b>           | <b>Q6VNS1</b>                                                                                                                                                                |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom


Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC


68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>MW</b>                  | 93 kDa                                                                                                                                  |
| <b>Tested applications</b> | IHC-P, WB                                                                                                                               |
| <b>Application notes</b>   | For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50                                                              |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                        |
| <b>Modifications</b>       | None                                                                                                                                    |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>NCBI</b>                | <b>Q6VNS1</b>                                                                                                                           |



Western blot analysis in mouse liver tissue lysates (35 ug/lane).



Mouse Ntrk3 Antibody immunohistochemistry analysis in formalin fixed and paraffin embedded mouse brain tissue followed by peroxidase conjugation of the secondary antibody and DAB staining.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558