

Product Datasheet

SERPINB4 Antibody (orb1270096)

Catalog Number	orb1270096
Category	Antibodies
Description	SERPINB4 Antibody
Target	SERPINB4
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Rabbit Ig
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	batch dependent
Buffer/Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Immunogen	This SPB4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 268-296 amino acids from the C-terminal region of human SPB4.
UniProt ID	P48594
MW	45 kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Tested applications	FC, WB
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50
Antibody Type	Primary Antibody
Modifications	None
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	P48594



Western blot analysis in K562 cell line lysates (35 ug/lane).



Flow cytometric analysis of K562 cells (right histogram) compared to a negative control cell (left histogram). FITC-conjugated goat-anti-rabbit secondary antibodies were used for the analysis.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558