

## Product Datasheet

### SARS-CoV-2 (COVID-19) ORF8 Antibody (orb1274101)

|                             |                                                                                                                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1274101                                                                                                                                                                                                                            |
| <b>Category</b>             | Antibodies                                                                                                                                                                                                                            |
| <b>Description</b>          | SARS-CoV-2 (COVID-19) ORF8 Antibody                                                                                                                                                                                                   |
| <b>Target</b>               | ORF8                                                                                                                                                                                                                                  |
| <b>Clonality</b>            | Polyclonal                                                                                                                                                                                                                            |
| <b>Conjugation</b>          | Unconjugated                                                                                                                                                                                                                          |
| <b>Reactivity</b>           | Virus                                                                                                                                                                                                                                 |
| <b>Form/Appearance</b>      | Liquid                                                                                                                                                                                                                                |
| <b>Concentration</b>        | batch dependent                                                                                                                                                                                                                       |
| <b>Buffer/Preservatives</b> | Antibody supplied in PBS.                                                                                                                                                                                                             |
| <b>Purification</b>         | Antibodies are affinity purified against the antigen and eluted with low pH.                                                                                                                                                          |
| <b>Immunogen</b>            | GST-SARS-CoV2 ORF8                                                                                                                                                                                                                    |
| <b>Tested applications</b>  | WB                                                                                                                                                                                                                                    |
| <b>Application notes</b>    | Detection Antibodies: 1 µg/ml Ab, 5% BSA in TBS-T, 4°C overnight Blocking: 5% BSA in TBS-T, 1h room temp Secondary Antibodies: 1:15,000 donkey anti-sheep, 1:15,000 donkey anti-rabbit, in TBS-T, 1h room temp Wash: TBS-T, 4x 10 min |
| <b>Antibody Type</b>        | Primary Antibody                                                                                                                                                                                                                      |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Storage**


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Dry Ice Shipping**

**Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.**

**Note**

For research use only



Antibody used at 1 ug/mL in WB.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558