

Product Datasheet

LAG-3 Single Domain Antibody [1A6] (orb1274187)


Description	LAG-3 Single Domain Antibody [1A6]
Species/Host	Camelus
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-P
Immunogen	LAG3 antibody was raised against a recombinant protein corresponding to amino acids 23-434 of human LAG3.
Target	LAG3
Preservatives	LAG-3 Antibody is supplied in PBS.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	LAG-3 antibody can be used for ELISA starting at 0.1 µg/mL. For immunohistochemistry starting at 1 µg/mL. For immunofluorescence start at 5 µg/mL. For flow cytometry at 1 µg/ml.
Isotype	sdAb
Clonality	Recombinant
Clone Number	1A6
Antibody Type	Primary Antibody

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Uniprot ID **P18627****NCBI** **NP_002277****Expiration Date** 12 months from date of receipt.

Immunohistochemistry of LAG-3 in human lymph node tissue with LAG-3 single domain antibody at 1 μ g/mL.



Flow cytometry analysis of LAG-3 overexpressing HEK293 cells using LAG-3 single domain antibody at 1 μ g/ml. Blue: Untransfected HEK293 cells. Yellow: LAG-3 overexpressing HEK293 cells.



Immunofluorescence of LAG-3 in transfected HEK293 cells with LAG-3 single domain antibody at 5 μ g/mL.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)