

Product Datasheet

PD-1 Single Domain Antibody [F12A8] (orb1274203)


Catalog Number	orb1274203
Description	PD-1 Single Domain Antibody [F12A8]
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	PD-1 antibody was raised against a recombinant protein corresponding to amino acids 25 - 167 of human PD-1.
Target	PDCD1
Preservatives	PD-1 Antibody is supplied in PBS.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	PD-1 antibody can be used for detection of PD-1 by ELISA and FACS (10 µg/ml).
Isotype	sdAb
Clonality	Recombinant
Clone Number	F1-2A8
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Uniprot ID **Q15116****NCBI** **NP_005009****Expiration Date** 12 months from date of receipt.

Titration ELISA of PD-1.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)