

Product Datasheet

Recombinant CD43 Antibody [SPN/1766R] (orb1274262)

Catalog Number	orb1274262
Description	Recombinant CD43 Antibody [SPN/1766R]
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Recombinant human protein was used as the immunogen for the recombinant CD43 antibody.
Target	SPN
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 1-2 ug/mlImmunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RTPrediluted IHC only format: incubate for 30 min at RT (1)Optimal dilution of the recombinant CD43 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.
Isotype	IgG, k


Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Clonality	Recombinant
Clone Number	SPN/1766R
Antibody Type	Primary Antibody
Uniprot ID	P16150
Expiration Date	12 months from date of receipt.



IHC testing of FFPE human tonsil with recombinant CD43 antibody (clone SPN/1766R). Required HIER: boil tissue sections in 10mM citrate buffer, pH 6, for 10-20 min followed by cooling at RT for 20 min.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)