

## Product Datasheet

### SLC9A7 Antibody (orb1277269)

|                             |                                                                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1277269                                                                                                                              |
| <b>Category</b>             | Antibodies                                                                                                                              |
| <b>Description</b>          | Rabbit polyclonal antibody to SLC9A7                                                                                                    |
| <b>Target</b>               | Disease                                                                                                                                 |
| <b>Clonality</b>            | Polyclonal                                                                                                                              |
| <b>Species/Host</b>         | Rabbit                                                                                                                                  |
| <b>Isotype</b>              | Rabbit IgG                                                                                                                              |
| <b>Conjugation</b>          | Unconjugated                                                                                                                            |
| <b>Reactivity</b>           | Human, Mouse                                                                                                                            |
| <b>Concentration</b>        | Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration.                                            |
| <b>Buffer/Preservatives</b> | Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.                                                                      |
| <b>Immunogen</b>            | Synthetic peptide of human SLC9A7                                                                                                       |
| <b>UniProt ID</b>           | <b>Q96T83</b>                                                                                                                           |
| <b>Tested applications</b>  | IHC                                                                                                                                     |
| <b>Dilution range</b>       | IHC: 1/50-1/200; ELISA: 1/2000-1/5000                                                                                                   |
| <b>Antibody Type</b>        | Primary Antibody                                                                                                                        |
| <b>Storage</b>              | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

|                        |                                 |
|------------------------|---------------------------------|
| <b>Note</b>            | For research use only           |
| <b>NCBI</b>            | <b>NP_115980.1</b>              |
| <b>Expiration Date</b> | 12 months from date of receipt. |

---

#### **Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### **Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)