

## Product Datasheet

### MRPL40 Antibody (orb1278320)

|                             |                                                                                              |
|-----------------------------|----------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1278320                                                                                   |
| <b>Category</b>             | Antibodies                                                                                   |
| <b>Description</b>          | Rabbit polyclonal antibody to MRPL40                                                         |
| <b>Target</b>               | MRPL40                                                                                       |
| <b>Clonality</b>            | Polyclonal                                                                                   |
| <b>Species/Host</b>         | Rabbit                                                                                       |
| <b>Isotype</b>              | Rabbit IgG                                                                                   |
| <b>Conjugation</b>          | Unconjugated                                                                                 |
| <b>Reactivity</b>           | Human                                                                                        |
| <b>Concentration</b>        | Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration. |
| <b>Buffer/Preservatives</b> | Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.                           |
| <b>Immunogen</b>            | Fusion protein of human MRPL40                                                               |
| <b>UniProt ID</b>           | <b>Q9NQ50</b>                                                                                |
| <b>MW</b>                   | 24 kDa                                                                                       |
| <b>Tested applications</b>  | IHC, WB                                                                                      |
| <b>Dilution range</b>       | WB: 1/1000-1/5000; IHC: 1/100-1/300; ELISA: 1/5000-1/10000                                   |
| <b>Antibody Type</b>        | Primary Antibody                                                                             |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)